

ANTERIOR VITRECTOMY SURGICAL EQUIPMENT

60-1500 Cuts Per Minute

GENERAL INFORMATION

Anterior Vitrectomy provides the following modes of operations

1. Suction
2. Cutting with Suction

The suction (Vacuum / Aspiration) can be varied from 0-500 mmHg and the cutting rate can be varied from 60 -1500 cuts/min by the incremental of 60 cuts/min.

FEATURES & SPECIFICATIONS

- Light weight and Portable
- Reusable collection bottle
- Easy to read 7 segment display
- Feather touch control keys for cut rate Up/Down
- Easy to use foot switch – cutting with suction or only suction
- 05-500 mmHg panel set venturi aspiration with digital display
- Pressure input: (5.5 – 6.5 Bar) or (80 – 95 PSI) - Air Compressor or Cylinder
- Electrical power: AC 100-240 Volts, 50/60 Hz, Single Phase
- Weight: 1.2 Kg / Dimensions: L 225mm x D 199 mm X H 72 mm



Manufactured by:



ICON MEDICARE

#304-A/2, TNHB, Phase VII, Brindavan Nagar,
HOSUR - 635 109. Tamilnadu, INDIA.

Phone : +91 4344 243313

E-mail : info@iconmedicare.com / www.iconmedicare.com

NEW DELHI
Mob.: +91 9944678009

MUMBAI
Mob.: +91 9323284971

BANGALORE
Mob.: +91 8220009934

CHENNAI
Mob.: +91 8144433313

KOLKATTA
Mob.: +91 9842191245

AHMEDABAD - Mob.: +91 9327065024

HYDERABAD - Mob.: +91 9000929890

ANTERIOR VITRECTOMY SURGICAL EQUIPMENT

RETIKARE[®]
FOR VISION

60-900 Cuts Per Minute

GENERAL INFORMATION

Anterior Vitrectomy provides the following modes of operations

1. Suction
2. Cutting with Suction

The suction (Vacuum / Aspiration) can be varied from 0-500 mmHg and the cutting rate can be varied from 60 -900 cuts/min by the incremental of 60 cuts/min.

FEATURES & SPECIFICATIONS

- Light weight and Portable
- Reusable collection bottle
- Easy to read 7 segment display
- Feather touch control keys for cut rate Up/Down
- Easy to use foot switch – cutting with suction or only suction
- 05-500 mmHg panel set venturi aspiration with digital display
- Pressure input: (5.5 – 6.5 Bar) or (80 – 95 PSI) - Air Compressor or Cylinder
- Electrical power: AC 100-240 Volts, 50/60 Hz, Single Phase
- Weight: 1.2 Kg / Dimensions: L 225mm x D 199 mm X H 72 mm



Manufactured by:



ICON MEDICARE

#304-A/2, TNHB, Phase VII, Brindavan Nagar,
HOSUR - 635 109. Tamilnadu, INDIA.

Phone : +91 4344 243313

E-mail : info@iconmedicare.com / www.iconmedicare.com

NEW DELHI
Mob.: +91 9944678009

MUMBAI
Mob.: +91 9323284971

BANGALORE
Mob.: +91 8220009934

CHENNAI
Mob.: +91 8144433313

KOLKATTA
Mob.: +91 9842191245

AHMEDABAD - Mob.: +91 9327065024

HYDERABAD - Mob.: +91 9000929890