

HIGH SPEED ANTERIOR VITRECTOMY SURGICAL EQUIPMENT

A-MAX (3000 Cuts Per Minute)



GENERAL INFORMATION

High Speed Anterior Vitrectomy System A-MAX provides the following modes of operations

VITRECTOMY

- Suction
- Cutting with Suction
- Reflux

WET-FIELD BIPOLAR COAGULATOR

Manufactured by:



INEST MEDICARE

#304-A/4, TNHB, Phase VII, Brindavan Nagar,
Hosur - 635 109. Tamil Nadu. INDIA

PH : +91 4344 - 243313 Email: info@inestmedicare.com

The Cutting rate can be varied from 60 -3000 cuts/min by the incremental of 60 cuts/min and the Suction (vacuum/aspiration) can be varied from 0-600 mmHg

FEATURES & SPECIFICATIONS

- 3000 Cuts/min
- 05-600 mmHg panel set venturi aspiration
- Reusable collection bottle with Fluid level Sensor
- Diathermy/Cautery & Forced Positive Reflux
- Light weight and Portable
- 20g,23g & 25g Compatible
- User Friendly 4.3" Touch TFT LCD
- Easy to use foot switch – cutting with suction, only suction & Reflux Function
- Pressure input: (6 – 7 Bar) or (90 – 105 PSI) - Air Compressor or Cylinder
- Weight: 2.5 Kg / Dimensions: L 225mm x D 199 mm X H 130 mm
- Electrical power: AC 100-240 Volts, 50/60 Hz, Single Phase



NEW DELHI
Mob.: +91 9944678009

MUMBAI
Mob.: +91 9323284971

BANGALORE
Mob.: +91 8220009934

CHENNAI
Mob.: +91 8144433313

KOLKATTA
Mob.: +91 9842191245

AHMEDABAD
Mob.: +91 9327065024

HYDERABAD
Mob.: +91 9000929890