

## LED COLD LIGHT SOURCE SINGLE

### RETIKARE - WHITE

#### GENERAL INFORMATION

Retikare-White LED Light Source, which offers the best bright pure white light for Ophthalmic surgical procedures. This light is offered with a first of its kind White Cold Light.

This LED light source is the most suitable alternative for halogen and xenon light sources.

#### FEATURES

- Long lifetime
- Compact and light weight
- PID controlled light intensity, ensuring longer life
- High brightness is comparable to xenon light source
- No harmful radiation of IR- infrared or UV-ultraviolet
- Quiet operation compared to conventional light sources

#### TECHNICAL SPECIFICATIONS

- Dimmer Control: Electronic
- Color Temperature: Natural Day Light
- Electrical power: AC 100-240 V, 50/60 Hz Single Phase
- Weight: 1.6 kg / Dimensions: L 225mm x D 199mm x H 72 mm



Manufactured by:



#### ICON MEDICARE

#304-A/2, TNHB, Phase VII, Brindavan Nagar,  
HOSUR - 635 109. Tamilnadu, INDIA.

Phone : +91 4344 243313

E-mail : [info@iconmedicare.com](mailto:info@iconmedicare.com) / [www.iconmedicare.com](http://www.iconmedicare.com)

**NEW DELHI**  
Mob.: +91 9944678009

**MUMBAI**  
Mob.: +91 9323284971

**BANGALORE**  
Mob.: +91 8220009934

**CHENNAI**  
Mob.: +91 8144433313

**KOLKATTA**  
Mob.: +91 9842191245

**AHMEDABAD - Mob.: +91 9327065024**

**HYDERABAD - Mob.: +91 9000929890**

## LED COLD LIGHT SOURCE DUAL

### RETIKARE - WHITE

#### GENERAL INFORMATION

Retikare-White LED Light Source, which offers the best bright pure white light for Ophthalmic surgical procedures. This light is offered with a first of its kind White Cold Light.

This LED light source is the most suitable alternative for halogen and xenon light sources.

#### FEATURES

- Long lifetime
- Compact and light weight
- PID controlled light intensity, ensuring longer life
- High brightness is comparable to xenon light source
- No harmful radiation of IR- infrared or UV-ultraviolet
- Quiet operation compared to conventional light sources
- Second LED source for stand-by / chandelier application

#### TECHNICAL SPECIFICATIONS

- Dimmer Control: Electronic
- Color Temperature: Natural Day Light
- Electrical power: AC 100-240 V, 50/60 Hz Single Phase
- Weight: 2.7 kg / L 290 mm x D 256 mm x H 62 mm



Manufactured by:



#### ICON MEDICARE

#304-A/2, TNHB, Phase VII, Brindavan Nagar,  
HOSUR - 635 109. Tamilnadu, INDIA.

Phone : +91 4344 243313

E-mail : [info@iconmedicare.com](mailto:info@iconmedicare.com) / [www.iconmedicare.com](http://www.iconmedicare.com)

**NEW DELHI**  
Mob.: +91 9944678009

**MUMBAI**  
Mob.: +91 9323284971

**BANGALORE**  
Mob.: +91 8220009934

**CHENNAI**  
Mob.: +91 8144433313

**KOLKATTA**  
Mob.: +91 9842191245

**AHMEDABAD - Mob.: +91 9327065024**

**HYDERABAD - Mob.: +91 9000929890**