

VIT PLUS 5000

VITRECTOMY CUT RATE ENHANCER WITH
SILICON OIL INJECTOR /REMOVER
& AIR INJECTOR (Optional)

The VIT PLUS 5000 allows enhanced High-Speed cutting (upto 5000 Cuts per Minute) with your host vitrectomy surgical system which has pneumatic cutter drive function.

VIT PLUS 5000 provides the following modes of operations.

- A. Vitrectomy Cutting (upto 5000CPM)
- B. Silicon Oil Injector / Remover
- C. Air injector (Optional)

FEATURES & SPECIFICATIONS

- Cutting Speed upto 5000 CPM
- User adjustable pressure setting (psi) to suit different types of cutters
- 4.3" TFT Touch Screen Display
- User Adjustable SOI – upto 80 psi / SOR - upto 600 mmHg
- SOI & SOR functions enabled through separate foot switch
- Optional Air Injector Function (Air/Fluid Exchanger)
- Light weight and Portable
- Pressure input: (7 – 8.5 Bar) or (100-125 psi) Compressed Air
- Electrical power: AC 100-240 Volts , 50/60 Hz, Single Phase
- Dimension: L 290 mm x D 256 mm x H 108 mm
- Weight:



Manufactured by:



ICON MEDICARE

#304-A/2, TNHB, Phase VII, Brindavan Nagar,
HOSUR - 635 109. Tamilnadu, INDIA.

Phone : +91 4344 243313

E-mail : info@iconmedicare.com / www.iconmedicare.com

NEW DELHI
Mob.: +91 9944678009

MUMBAI
Mob.: +91 9323284971

BANGALORE
Mob.: +91 8220009934

CHENNAI
Mob.: +91 8144433313

KOLKATTA
Mob.: +91 9842191245

AHMEDABAD - Mob.: +91 9327065024

HYDERABAD - Mob.: +91 9000929890

RETIKARE I/P

AIR INJECTOR / GAS FORCED INFUSION

The RETIKARE I/P unit allows you to use as Air Injector Function or Infusion Pressure (Air Forced Infusion) with Temponade for both Anterior/Posterior Surgical Procedures

RETIKARE I/P provides you the following modes of operations.

- A. Air Injector Function (or)
- B. Air Forced Infusion Pressure with Temponade Function

FEATURES & SPECIFICATIONS

- Positive Air Pressure Output upto 100mmHg
- Inbuild Air Pump (Eliminates External Air Pressure source)
- Temponade Function enable/disable through separate Foot switch
- Used for both Anterior/Posterior Surgical Procedures
- Easy to read 7 segment display
- Feather touch control keys for cut rate Up/Down
- Light weight and Portable
- Electrical power: AC 100-240 Volts, 50/60 Hz, Single Phase
- Weight: 1.2 Kg
- Dimensions: L 225mm x D 199 mm X H 72 mm



Manufactured by:



ICON MEDICARE

#304-A/2, TNHB, Phase VII, Brindavan Nagar,
HOSUR - 635 109. Tamilnadu, INDIA.

Phone : +91 4344 243313

E-mail : info@iconmedicare.com / www.iconmedicare.com

NEW DELHI
Mob.: +91 9944678009

MUMBAI
Mob.: +91 9323284971

BANGALORE
Mob.: +91 8220009934

CHENNAI
Mob.: +91 8144433313

KOLKATTA
Mob.: +91 9842191245

AHMEDABAD - Mob.: +91 9327065024

HYDERABAD - Mob.: +91 9000929890